

Production of Bioactive Isothiocyanates by Bacterial Myrosinase

Juan Isaiz,¹ Xiaoji Liu,² and Michael J. Miller²

Moraine Valley Community College, Palos Hills, Illinois¹

Department of Food Science and Human Nutrition, University of Illinois at Urbana-Champaign²

PRECS Phenotypic Plasticity Research Experience for Community College Students

Introduction

- The focus of the project is on the hydrolysis of glucosinolates (GSLs) into isothiocyanates (ITC) by gut microbial myrosinase.
- ITCs are bioactive products found in several plants belonging to the *Brassicaceae* family, and are produced upon disruption of the tissue by herbivores or pathogens.
- ITCs are normally formed within plants as a defense against herbivores but in humans they have shown anti-cancer and anti-inflammatory properties.¹




Figure 1: Glucosinolate hydrolysis and its various products.³

Background

- Although we know gut microbes can hydrolyze GSL into ITC, the mechanism is not known.
- Previous studies have shown that a probiotic strain of bacteria *Escherichia coli* Nissle 1917 is able to hydrolyze GSL without any detectable ITCs.⁵
- *E. coli* Nissle 1917, isolated from a German soldier during World War I, is now being used to treat gastrointestinal disorders.⁴
- *Citrobacter* Wye1 was isolated from a UK soil by enrichment using sinigrin.
- With improved detection methods, the aim of this experiment is to explore whether the bacterial myrosinase β -glucosidase found within the genome of *E. coli* Nissle 1917 can hydrolyze GSL into detectable ITCs.

Question


- Will the microbial myrosinase hydrolyze the GSL into bioactive ITC?

Hypothesis

- Microbial myrosinase is expected to produce bioactive ITCs but not as much as plant myrosinase.

Methods

- For this experiment we have chosen two candidate myrosinases from the gut microbe *E. coli* Nissle 1917 and the soil microbe *Citrobacter*. *Citrobacter* is the positive control in this experiment.
- We have the genes from these two candidate myrosinases synthesized (ThermoFisher) and cloned into expression plasmids pRSET A.



- We transformed *E. coli* BL21(DE3) with the plasmids. We plated them on LB inoculated with Ampicillin, positive clones were picked and checked by DNA electrophoresis.
- The overexpression of microbial myrosinases was induced by IPTG. The myrosinases were purified by Ni-NTA column.
- The hydrolysis of GSL, sinigrin, was examined by cyclocondensation-HPLC.
- The hydrolysis of 4-nitrophenyl β -d-glucopyranoside (4NPG) was examined by A405nm

Results

Fig 1. DNA Gel Electrophoresis confirming the sizes of plasmids cloned into *E. coli* BL21(DE3)




Fig 1. Transformation of plasmids with genes encoding myrosinase from *Citrobacter* sp. Wye1 myrosinase and 6-phospho-beta-glucosidase gene from *E. coli* Nissle 1917 into *E. coli* BL21(DE3).

- The overexpression of microbial myrosinases was induced by IPTG. The myrosinases were purified by Ni-NTA column.

Fig 2. Purification of Nissle & *Citrobacter* myrosinase by Ni-NTA columns and Amicon® Ultra Centrifugal Filters (30k MWCO)




Fig 2. Purified Citro_WYE1 and Nissle_15645 (with 6x His-tag) visualized on 12% SDS-PAGE gel.

Results con't

- The hydrolysis of GSL, sinigrin, was examined by cyclocondensation-HPLC.

Fig 3. HPLC chromatogram showing the ITC production




Fig 3. HPLC chromatogram showing the ITC production by *E. coli* BL21(DE3) transformed with pRSET_A plasmid expressing *Citrobacter* WYE1 myrosinase⁴²

Fig 4. Sinigrin hydrolysis by purified microbial myrosinase candidates




Fig 4. The myrosinase from *Citrobacter* produced detectable ITCs but not the myrosinase from *E. coli* Nissle.

- The hydrolysis of 4-nitrophenyl β -d-glucopyranoside (4NPG) was examined by A405nm

Fig 5. 4NPG hydrolysis by purified microbial myrosinase candidates




Fig 5. The myrosinase from *Citrobacter* produced detectable p-nitrophenol over time but not the myrosinase from *E. coli* Nissle.

Conclusions

- We have first-documented, direct evidence that a microbial myrosinase can produce detectable ITCs from GSLs.
- It is possible that the Nissle myrosinase is still a functioning enzyme, but only capable of hydrolyzing phosphorylated substrate, not GSLs/4NPG.

Future Work

- Use knock-out gene library strategy to identify the gene(s) encoding enzyme(s) that are responsible for the production of ITCs.

References

1. Angelino, Donato et al. "Myrosinase-Dependent and -independent Formation and Control of Isothiocyanate Products of Glucosinolate Hydrolysis." *Frontiers in Plant Science* 6 (2015): 831. *PMC*. Web. 9 July 2018.
2. Albaser, A.; Kazana, E.; Bennett, M. H.; Cebeci, F.; Luang-In, V.; Spanu, P. D.; Rossiter, J. T. Discovery of a bacterial glycoside hydrolase family 3 (GH3) β -glucosidase with myrosinase activity from a *Citrobacter* strain isolated from soil. *Journal of Agricultural and Food Chemistry* 2016, 64, 1520–1527, doi:10.1021/acs.jafc.5b05381.
3. F. Vaughn, Steven & Berhow, Mark. (2005). Glucosinolate hydrolysis products from various plant sources: pH effects, isolation, and purification. *Industrial Crops and Products*. 21, 193-202. 10.1016/j.indcrop.2004.03.004.
4. Liu, Xiaoji et al. "Dietary Broccoli Alters Rat Cecal Microbiota to Improve Glucoraphanin Hydrolysis to Bioactive Isothiocyanates." *Nutrients* 9.3 (2017): 262. *PMC*. Web. 9 July 2018.
5. Ulrich Sonnenborn; *Escherichia coli* strain Nissle 1917—from bench to bedside and back: history of a special *Escherichia coli* strain with probiotic properties, *FEMS Microbiology Letters*, Volume 363, Issue 19, 1 October 2016, fnw212, <https://doi.org/10.1093/femsle/fnw212>
6. Mullaney J. A., Kelly W. J., McGhie T. K., Ansell J., Heyes J. A. (2013). Lactic acid bacteria convert glucosinolates to nitriles efficiently yet differently from enterobacteriaceae. *J Agric. Food Chem.* 61, 3039–3046. 10.1021/jf305442j

Acknowledgments

Financial support was provided by the National Science Foundation under grant #NSF REU 1559908/1559929, as part of the Phenotypic Plasticity Research Experience for Community College Students, through the University of Illinois at Urbana-Champaign Institute for Genomic Biology and Parkland College. <http://precs.igb.illinois.edu/>

I would like to thank Luis Ibarra, Max Holle, Shannon Rezac, Dr. Liu and Dr. Miller for volunteering their time this summer to mentor me. Thanks also to everyone at the Miller for their kindness and support. A huge thanks to the projects PIs Dr. Nathan Schroeder and Dr. C. Britt Carlson for creating the program. Thanks also to the technical and support staff at the Institute for Genomic Biology.